

**SOP:** Propagation of HepG2  
**Date modified:** 9/5/2008  
**Modified by:** J. Goldy/M. Dorschner

### **Ordering Information**

HepG2 can be ordered from ATCC as a frozen ampoule.

Name: HepG2, hepatocellular carcinoma  
ATCC #: HB-8065

### **Notes:**

This is an adherent cell line.

### **Materials List**

1. MEM with 2mM L-glutamine (Cellgro Cat# 10-010-CM)
2. Fetal Bovine Serum (Cellgro Cat# 35-016-CV)
3. Sodium Bicarbonate (Cellgro Cat# 25-035-CI)
4. Sodium pyruvate (Cellgro Cat# 25-000-CI)
5. Non-essential amino acids (Cellgro Cat# 25-025-CI)
6. T75 & T225 culture flasks
7. Graduated pipets (1, 5, 25mL)
8. Penicillin-Streptomycin Solution (100X) (Cellgro Cat# 30-002-CI)
9. Hemocytometer
10. Micropipet w/ P20 tips
11. Microscope

### **Growth Medium for HepG2**

MEM with 2mM L-glutamine  
Non-essential amino acids  
Sodium Pyruvate 1mM  
10% FBS  
Pen-Strep (1X)  
Sodium Bicarbonate 1.5g/L

### **Procedure**

#### **A. Receipt of frozen cells and starting cell cultures.**

- 1) Immediately place frozen cells in liquid nitrogen storage incubator.
- 2) Quickly thaw ampoule in 37°C water bath
- 3) Transfer thawed cells to a T75 flask with 40mLs of warm growth media.
- 4) Allow cells to recover over night in 37°C, 5% CO<sub>2</sub> humidified incubator.
- 5) Pour off medium the next day, replace with fresh medium and return to incubator.

#### **B. Sub-culture**

- 1) Propagate cells until density reaches 70-80% confluence.
- 2) Decant medium.
- 3) Wash cells with warm 1X PBS.
- 4) Add 8mLs of Accutase and return to incubator for 10-15 minutes.

- 5) Immediately remove cells and pellet at 500 xg for 5 minutes (4°C).
- 6) Wash cells 2X with 1X PBS.
- 7) Gently re-suspend cell pellet in warm medium.
- 8) Perform 1:4 to 1:6 cell split as needed.
- 9) Record each subculture event as a passage.

**C. Maintenance**

- 1) Change media the day after seeding and every 2-3 days thereafter.  
Use ~50mLs of medium per T225 flasks.

**D. Harvest**

- 1) Do not use cells that have been passed more than 8 times.
- 2) Remove cells from flasks according to protocol described above under 'subculturing'.
- 3) Examine viability using trypan blue staining (SOP).