

SOP: Propagation of MCF7
Date modified: 12/14/09
Modified by: T. Canfield (UW)

Ordering Information

MCF7 can be ordered from ATCC as a frozen ampoule.

Name: MCF7, mammary gland adenocarcinoma
ATCC #: HTB-22

Notes:

This is an adherent cell line.

Materials List

1. Eagle's MEM with Earle's salts and L-Glutamine (Cellgro Cat# 10-010-CM)
2. Characterized Fetal Bovine Serum (HyClone Cat# SH30071)
3. Non-essential Amino Acids, 100X (Invitrogen, Cat#11140-050)
4. Insulin from bovine pancreas, cell-culture tested, 200X (Sigma-Aldrich Cat# I6634, prepared in acidified water per manufacturer's specifications)
5. T75 and T225 culture flasks
6. Graduated pipets (1, 5, 10, 25, 50mL)
7. Penicillin-Streptomycin Solution, 200X (Cellgro, Cat# 30-001-CI)
8. Phosphate Buffered Saline (1X PBS) (prepared from 10X stock Cellgro, Cat# 46-013-CM by dilution with sterile deionized water)
9. Accutase Enzyme Cell Detachment Medium (EBioscience, Cat# 00-455)
10. Freezing medium (growth medium containing 6% DMSO)
11. DMSO, Hybri-Max (Sigma-Aldrich Cat# D2650)
12. Cryovials (Nunc Cat# 368632)
13. Hemocytometer
14. Micropipet w/ P20 tips
15. Microscope

Growth Medium for MCF7

Eagle's MEM with Earle's salts and L-Glutamine
Non-essential Amino Acids (1X)
0.01 mg/ml Bovine Insulin
10% FBS
Pen-Strep (1X)

Procedure

A. Receipt of Frozen Cells and Starting Cell Culture

- 1) Immediately place frozen cells in liquid nitrogen freezer storage until ready to culture.
- 2) Quickly thaw ampoule in 37°C water bath.
- 3) Transfer thawed cells to a T75 flask with 20ml of warm growth media.
- 4) Allow cells to recover overnight in 37°C, 5% CO₂ humidified incubator.
- 5) The next morning, the diluted DMSO-containing shipping/cryopreservation medium is aspirated from the cell layer and replaced with fresh medium.

B. Sub-culture

- 1) Propagate cells until density reaches 70-80% confluence.
- 2) Aspirate medium.
- 3) Wash cells with warm 1X PBS.
- 4) Add 15mLs of Accutase and return to incubator for 10-15 minutes or until cells detach.
- 5) Immediately remove cells, rinse flask with warm 1X PBS to collect residual cells, and pellet at 500 X g for 5 minutes (4°C).
- 6) Gently re-suspend cell pellet in warm medium.
- 7) Perform 1:3 to 1:6 cell split as needed.
- 8) Record each subculture event as a passage.

C. Maintenance and Generation of Seed Stocks

- 1) Change medium the day after seeding and every 2-3 days thereafter. Use ~50mLs of medium per T225 flask.
- 2) Following first or second passage after receipt of cells and with sufficient number of cells to continue maintenance and expansion, a portion of the flasks should be sub-cultured using Accutase as above under “Sub-culture” and the cell pellet resuspended in freezing medium.
- 3) Cells are dispensed into cryovials (2 million cells per 1 mL aliquot) and frozen in a -80°C isopropanol cryo-freezing container overnight.
- 4) Cryovials are transferred the next day to liquid nitrogen freezer for long-term storage.

D. Harvest

- 1) Passage cells until the desired number of cells is reached.
- 2) Pellet cells and rinse with 1X PBS as described above in “Sub-culture and Maintenance”.
- 3) Examine viability using Trypan blue staining (SOP TP-7).