

**SOP:** Propagation of Mouse A20 Lymphoblasts  
**Date modified:** 01/18/2011  
**Modified by:** T. Canfield (UW)

### **Ordering Information**

Mouse A20 Lymphoblasts can be ordered from ATCC as a frozen ampoule. This is a mouse suspension lymphoblast cell line derived from a reticulum cell sarcoma (B cell lymphoma).

Name: A20—Mouse Reticulum Cell Sarcoma Lymphoblast  
ATCC #: TIB-208

### **Materials List**

1. RPMI 1640 with 2mM L-glutamine adjusted to contain 1.5g/L sodium bicarbonate, 4.5g/L glucose, 10mM HEPES, and 1.0mM sodium pyruvate Medium (ATCC, Cat# 30-2001)
2. Characterized Fetal Bovine Serum (HyClone, Cat# SH30071)
3. T75, T225 tissue culture flasks
4. Corning conical centrifuge tubes (15mL and 50mL)
5. Graduated pipets (1, 5, 10, 25, 50mL)
6. Penicillin-Streptomycin Solution (200X) (Cellgro, Cat# 30-001-CI)
7. Phosphate Buffered Saline (1X PBS) (Cellgro, Cat# 21-040-CM)
8. 2-Mercaptoethanol, liquid, cell culture tested (Sigma-Aldrich, Cat# M7522)
9. Freezing Medium (growth medium containing 5% DMSO)
10. DMSO, Hybri-Max (Sigma-Aldrich, Cat# D2650)
11. Cryovials (Nunc, Cat# 368632)
12. Cryo 1°C Freezing Container (Nalgene, Cat# 5100-0001)
13. Eppendorf Centrifuge 5810R
14. Revco UltimaII -80°C Freezer
15. Thermolyne Locator 4 Liquid Nitrogen Freezer
16. Hemocytometer
17. Micropipet w/ P20 tips
18. Microscope

### **Growth Medium for Mouse A20 Lymphoblasts**

RPMI 1640 with 2mM L-glutamine adjusted to contain 1.5g/L sodium bicarbonate, 4.5g/L glucose, 10mM HEPES, and 1.0mM sodium pyruvate Medium  
10% Characterized Fetal Bovine Serum  
0.05mM 2-Mercaptoethanol  
Pen-Strep (1X)

### **Procedure**

#### **A. Receipt of Frozen Cells and Starting Cell Culture**

- 1) Immediately place frozen cells in liquid nitrogen storage until ready to culture.
- 2) When ready to start cell culture, quickly thaw ampoule in a 37°C water bath.
- 3) As soon as ice crystals disappear, swab outside surface of the ampoule with 70% ethanol, and then transfer the vial contents to a 15mL conical centrifuge tube containing 9mL of complete growth medium.
- 4) Centrifuge at 125 x g (4°C) for 7 minutes.
- 5) Resuspend the cell pellet in complete pre-warmed growth medium at a cell concentration of  $2 \times 10^5$  cells/mL in a T75 flask and grow in a 37°C, 5% CO<sub>2</sub> humidified incubator.

**B. Sub-culture and Maintenance**

1. Take cell counts with a hemocytometer every 24 hours to maintain the culture at a cell density between  $1 \times 10^5$  cells/mL and  $1 \times 10^6$  cells/mL. The cells have a doubling time of 18 hours and **the concentration of cells should not exceed  $1 \times 10^6$  cells/mL.**
2. Add fresh warm medium when appropriate to maintain cell density and expand the culture to the desired number of cells.
3. Record each subculture event as a passage.

**C. Generation of Seed Stocks**

1. At an early stage of expansion and with sufficient number of cells to continue maintenance, a small portion of the cells should be set aside as a seed stock, if needed.
2. Amount of cells for the seed stock should be placed in a conical centrifuge tube and centrifuged at  $500 \times g$  ( $4^\circ\text{C}$ ) for 5 minutes.
3. Aspirate supernatant and resuspend the cell pellet in 1X PBS to wash. Centrifuge again under same conditions.
4. Resuspend the cell pellet in freezing medium (growth medium containing 5% DMSO) at a concentration yielding 2 million cells per 1mL aliquot.
5. Dispense 1mL cell suspension per cryovial. Place cryovials in a Nalgene Cryo  $1^\circ\text{C}$  freezing container and store overnight at  $-80^\circ\text{C}$ .
6. Cryovials are transferred the next day to liquid nitrogen freezer for long-term storage.

**D. Harvest**

1. Passage cells until the desired number of cells for experimentation is reached in a logarithmic growth phase.
2. Pellet cells and rinse with 1X PBS as in “Generation of Seed Stocks” section.
3. Examine viability using Trypan blue staining (SOP TP-7).